

Scalpels



comknife®

Dura Blade
Straight Blade
Arachnoid Blade
Curved Blade
L-Hook Blade
Disposable Microsurgical Scalpels



Contents

comknife® Series

Disposable Microsurgical Scalpels (Neurosurgical Knife)	...04
Disposable Microsurgical Scalpels (Vessel Knife)	...06
Disposable Microsurgical Scalpels (Ophthalmic Knife)	...10

Neurosurgical Knife

Blade



Dura Blade

DT-Z-352

This blade is specially designed for dura mater. Its curved and blunt guard protects the tissues underneath during incision of dura mater or spinal dura mater.



Dura Blade

DT-W-353

This blade is specially designed with bend angle for dura mater. Its curved and blunt guard protects the tissues underneath during incision of dura mater or spinal dura mater.



Arachnoid Blade

DT-Z-354

This Blade is designed for arachnoid blade. The 0.3mm depth-limit design on the edge of circular blade prevents damage to deeper tissues.

The 2.0mm wide edge reveals excellent view of the incision.



Straight Blade

DT-Z-356

The ultra sharp 25° straight beveled blade is designed with a pointed tip at the front end. It is suitable for puncturing and cutting of micro-incision in narrow parts.



Curved Blade

DT-Z-357

The ultra sharp curved blade is with unique arc design can be used to hook tissue for cutting.



L-Hook Blade

DT-Z-358

The unique L-hook shaped blade with a small sickle design can be used to hook tissue for sharp piercing, and sharp cutting.

Model	Product Name	Material	Blade Length	Weight (Unit)	Secondary Package	Shipping Package
DT-Z-352	Dura Blade	Stainless Steel (30Cr13) + ABS + Titanium alloy (TC4)	21 mm	0.503 g	5 Pcs./Box	300 Pcs. / Ctn.
DT-W-353	Dura Blade	Stainless Steel (30Cr13) + ABS + Titanium alloy (TC4)	20 mm	0.503 g	5 Pcs./Box	300 Pcs. / Ctn.
DT-Z-354	Arachnoid Blade	Stainless Steel (30Cr13) + ABS + Titanium alloy (TC4)	10 mm	0.503 g	5 Pcs./Box	300 Pcs. / Ctn.
DT-Z-356	Straight Blade	Stainless Steel (30Cr13) + ABS + Titanium alloy (TC4)	18 mm	0.392 g	5 Pcs./Box	300 Pcs. / Ctn.
DT-Z-357	Curved Blade	Stainless Steel (30Cr13) + ABS + Titanium alloy (TC4)	18 mm	0.387 g	5 Pcs./Box	300 Pcs. / Ctn.
DT-Z-358	L-Hook Blade	Stainless Steel (30Cr13) + ABS + Titanium alloy (TC4)	18 mm	0.395 g	5 Pcs./Box	300 Pcs. / Ctn.

Curved Handle



Model	Material	Length	Weight (Unit)	Secondary Package	Shipping Package
DB-W-361	Titanium Alloy (TC4)	185 mm	15 g	1 Piece / Box	20 Pcs. / Ctn.
DB-W-362	Titanium Alloy (TC4)	210 mm	15.368 g	1 Piece / Box	20 Pcs. / Ctn.
DB-W-363	Titanium Alloy (TC4)	310 mm	15 g	1 Piece / Box	20 Pcs. / Ctn.

Straight Handle



Model	Material	Length	Weight (Unit)	Secondary Package	Shipping Package
DB-Z-364	Titanium Alloy (TC4)	100 mm	13 g	1 Piece / Box	20 Pcs. / Ctn.
DB-Z-365	Titanium Alloy (TC4)	130 mm	13 g	1 Piece / Box	20 Pcs. / Ctn.
DB-Z-366	Titanium Alloy (TC4)	150 mm	13 g	1 Piece / Box	20 Pcs. / Ctn.
DB-Z-367	Titanium Alloy (TC4)	175 mm	13.876 g	1 Piece / Box	20 Pcs. / Ctn.

Puncture Type Handle



Model	Material	Length	Weight (Unit)	Secondary Package	Shipping Package
DB-W-368	Titanium Alloy (TC4)	100 mm	/	1 Piece / Box	10 Pcs. / Ctn.
DB-W-369	Titanium Alloy (TC4)	200 mm	/	1 Piece / Box	10 Pcs. / Ctn.

Vessel Knife

Integrated Type



This series is specially designed for fine dissection and dissection of blood vessels and soft tissues in microsurgery, which not only ensures safety, but also makes operation much easier.

Stripping Knife

- Double edge design
- Sharp and convenient
- Swan neck design to reduce the blocking of vision under microscope during the operation.



Model	Material		Blade Length	Handle Length	Weight (Unit)	Secondary Package	Shipping Package
	Blade	Handle					
BL-B-117	Stainless Steel (30Cr13)	ABS	12.5 mm	176 mm	5.618 g	1 Piece / Box	100 Pcs. / Ctn.
BL-B-118	Stainless Steel (30Cr13)	ABS	14.5 mm	176 mm	5.618 g	1 Piece / Box	100 Pcs. / Ctn.
BL-B-119	Stainless Steel (30Cr13)	ABS	16 mm	176 mm	5.618 g	1 Piece / Box	100 Pcs. / Ctn.

Dissection Knife



Model	Material		Blade Length	Handle Length	Weight (Unit)	Secondary Package	Shipping Package
	Blade	Handle					
PQ-L-111	Stainless Steel (30Cr13)	ABS	3 mm	186 mm	5.408 g	1 Piece / Box	100 Pcs. / Ctn.
PQ-L-112	Stainless Steel (30Cr13)	ABS	4.5 mm	186 mm	5.408 g	1 Piece / Box	100 Pcs. / Ctn.
PQ-L-113	Stainless Steel (30Cr13)	ABS	6 mm	186 mm	5.408 g	1 Piece / Box	100 Pcs. / Ctn.



Model	Material		Blade Length	Angle	Weight (Unit)	Secondary Package	Shipping Package
	Blade	Handle					
PQ-J-114	Stainless Steel (30Cr13)	ABS	3 mm	35°	5.544 g	1 Piece / Box	100 Pcs. / Ctn.
PQ-J-115	Stainless Steel (30Cr13)	ABS	4.5 mm	35°	5.544 g	1 Piece / Box	100 Pcs. / Ctn.
PQ-J-116	Stainless Steel (30Cr13)	ABS	6 mm	35°	5.544 g	1 Piece / Box	100 Pcs. / Ctn.

Vessel Knife

Separated Type



Handle



Model	Material	Length	Angle	Weight (Unit)	Secondary Package	Shipping Package
DB-Z	304 Stainless Steel	177 mm	0°	26.760 g	1 Piece / Box	20 Pcs. / Ctn.
DB-W	304 Stainless Steel	175 mm	10°	26.446 g	1 Piece / Box	20 Pcs. / Ctn.

Stripping Knife



Model	Material	Blade Length	Weight (Unit)	Secondary Package	Shipping Package
BL-B-1	Stainless Steel (30Cr13) + ABS	12.5 mm	0.487 g	5 Pcs. / Box	300 Pcs. / Ctn.
BL-B-2	Stainless Steel (30Cr13) + ABS	14.5 mm	0.487 g	5 Pcs. / Box	300 Pcs. / Ctn.
BL-B-3	Stainless Steel (30Cr13) + ABS	16 mm	0.487 g	5 Pcs. / Box	300 Pcs. / Ctn.

Dissection Knife



0°

0° scalpel + 0° handle



-10°

0° scalpel - 10° handle



10°

0° scalpel + 10° handle



Model	Material	Blade Length	Angle	Weight (Unit)	Secondary Package	Shipping Package
PQ-L-1	Stainless Steel (30Cr13) + ABS	3 mm	0°	0.186 g	5 Pcs. / Box	300 Pcs. / Ctn.
PQ-L-2	Stainless Steel (30Cr13) + ABS	4.5 mm	0°	0.186 g	5 Pcs. / Box	300 Pcs. / Ctn.
PQ-L-3	Stainless Steel (30Cr13) + ABS	6 mm	0°	0.186 g	5 Pcs. / Box	300 Pcs. / Ctn.



25°

35° scalpel - 10° handle



35°

35° scalpel + 0° handle



45°

35° scalpel + 10° handle



Model	Material	Blade Length	Angle	Weight (Unit)	Secondary Package	Shipping Package
PQ-J-1	Stainless Steel (30Cr13) + ABS	3 mm	35°	0.186 g	5 Pcs. / Box	300 Pcs. / Ctn.
PQ-J-2	Stainless Steel (30Cr13) + ABS	4.5 mm	35°	0.186 g	5 Pcs. / Box	300 Pcs. / Ctn.
PQ-J-3	Stainless Steel (30Cr13) + ABS	6 mm	35°	0.186 g	5 Pcs. / Box	300 Pcs. / Ctn.

Ophthalmic Knife

Slit Knife


Usage: to make transparent corneal incision in cataract surgery



Model	Material		Bend Angle	Edge Width	Blade Thickness	Handle Length	Weight (Unit)	Secondary Package
	Blade	Handle						
CC-L-220	Stainless Steel (30Cr13)	ABS	Straight	1.0 mm	0.15 mm	117 mm	2.972 g	6 Pcs. / Box
CC-L-222	Stainless Steel (30Cr13)	ABS	Slit Angled 45°	1.0 mm	0.15 mm	117 mm	2.972 g	6 Pcs. / Box
CC-L-223	Stainless Steel (30Cr13)	ABS	Slit Angled 45°	1.1 mm	0.15 mm	117 mm	2.972 g	6 Pcs. / Box
CC-L-224	Stainless Steel (30Cr13)	ABS	Slit Angled 45°	1.2 mm	0.15 mm	117 mm	2.972 g	6 Pcs. / Box
CC-L-225	Stainless Steel (30Cr13)	ABS	Slit Angled 45°	1.3 mm	0.15 mm	117 mm	2.972 g	6 Pcs. / Box
CC-L-226	Stainless Steel (30Cr13)	ABS	Slit Angled 45°	1.4 mm	0.15 mm	117 mm	2.972 g	6 Pcs. / Box
CC-L-227	Stainless Steel (30Cr13)	ABS	Slit Angled 45°	1.5 mm	0.15 mm	117 mm	2.972 g	6 Pcs. / Box
CC-L-228	Stainless Steel (30Cr13)	ABS	Slit Angled 45°	3.0 mm	0.15 mm	117 mm	2.972 g	6 Pcs. / Box
CC-L-229	Stainless Steel (30Cr13)	ABS	Slit Angled 45°	1.6 mm	0.15 mm	117 mm	2.972 g	6 Pcs. / Box
CC-L-230	Stainless Steel (30Cr13)	ABS	Slit Angled 45°	1.7 mm	0.15 mm	117 mm	2.972 g	6 Pcs. / Box
CC-L-231	Stainless Steel (30Cr13)	ABS	Slit Angled 45°	1.8 mm	0.15 mm	117 mm	2.972 g	6 Pcs. / Box
CC-L-232	Stainless Steel (30Cr13)	ABS	Slit Angled 45°	1.9 mm	0.15 mm	117 mm	2.972 g	6 Pcs. / Box
CC-L-233	Stainless Steel (30Cr13)	ABS	Slit Angled 45°	2.0 mm	0.15 mm	117 mm	2.972 g	6 Pcs. / Box
CC-L-234	Stainless Steel (30Cr13)	ABS	Slit Angled 45°	2.1 mm	0.15 mm	117 mm	2.972 g	6 Pcs. / Box
CC-L-235	Stainless Steel (30Cr13)	ABS	Slit Angled 45°	2.2 mm	0.15 mm	117 mm	2.972 g	6 Pcs. / Box
CC-L-236	Stainless Steel (30Cr13)	ABS	Slit Angled 45°	2.3 mm	0.15 mm	117 mm	2.972 g	6 Pcs. / Box
CC-L-237	Stainless Steel (30Cr13)	ABS	Slit Angled 45°	2.4 mm	0.15 mm	117 mm	2.972 g	6 Pcs. / Box
CC-L-238	Stainless Steel (30Cr13)	ABS	Slit Angled 45°	2.5 mm	0.15 mm	117 mm	2.972 g	6 Pcs. / Box
CC-L-239	Stainless Steel (30Cr13)	ABS	Slit Angled 45°	2.65 mm	0.15 mm	117 mm	2.972 g	6 Pcs. / Box
CC-L-240	Stainless Steel (30Cr13)	ABS	Slit Angled 45°	2.75 mm	0.15 mm	117 mm	2.972 g	6 Pcs. / Box
CC-L-241	Stainless Steel (30Cr13)	ABS	Slit Angled 45°	2.85 mm	0.15 mm	117 mm	2.972 g	6 Pcs. / Box
CC-L-242	Stainless Steel (30Cr13)	ABS	Slit Angled 45°	2.95 mm	0.15 mm	117 mm	2.972 g	6 Pcs. / Box
CC-L-243	Stainless Steel (30Cr13)	ABS	Slit Angled 45°	3.2 mm	0.15 mm	117 mm	2.972 g	6 Pcs. / Box
CC-L-244	Stainless Steel (30Cr13)	ABS	Slit Angled 45°	3.5 mm	0.15 mm	117 mm	2.972 g	6 Pcs. / Box

Stab Knife

Usage: lateral corneal incision and scleral tunnel incision in cataract surgery



Model	Material		Bend Angle	Edge Width	Blade Thickness	Handle Length	Weight (Unit)	Secondary Package
	Blade	Handle						
CQ-L-200	Stainless Steel (30Cr13)	ABS	Double Blade, 0°	11 mm	0.15 mm	117 mm	2.940 g	6 Pcs. / Box
CQ-L-201	Stainless Steel (30Cr13)	ABS	Double Blade, 22.5°	11 mm	0.15 mm	117 mm	2.940 g	6 Pcs. / Box
CQ-L-202	Stainless Steel (30Cr13)	ABS	Double Blade, 15°	11 mm	0.15 mm	117 mm	2.940 g	6 Pcs. / Box
CQ-L-203	Stainless Steel (30Cr13)	ABS	Double Blade, 30°	11 mm	0.15 mm	117 mm	2.940 g	6 Pcs. / Box
CQ-L-204	Stainless Steel (30Cr13)	ABS	Double Blade, 45°	11 mm	0.15 mm	117 mm	2.940 g	6 Pcs. / Box

Crescent Knife

Usage: scleral tunnel incision for trabeculectomy in cataract surgery



Model	Material		Bend Angle	Edge Width	Blade Thickness	Handle Length	Weight (Unit)	Secondary Package
	Blade	Handle						
SD-L-210	Stainless Steel (30Cr13)	ABS	Crescent Bevel Up, 45°	1.6 mm	0.15 mm	117 mm	2.951 g	6 Pcs. / Box
SD-L-211	Stainless Steel (30Cr13)	ABS	Crescent Bevel Up, 45°	2.2 mm	0.15 mm	117 mm	2.951 g	6 Pcs. / Box
SD-L-212	Stainless Steel (30Cr13)	ABS	Crescent Bevel Up, 45°	2.6 mm	0.15 mm	117 mm	2.951 g	6 Pcs. / Box
SD-L-213	Stainless Steel (30Cr13)	ABS	Crescent Bevel Up, 60°	2.6 mm	0.15 mm	117 mm	2.951 g	6 Pcs. / Box

MVR Knife

Usage: Scleral incision, side port creation
 Corneal incision in cataract surgery
 Minor and precise incision in vitreous surgery and vitrectomy of
 posterior segment of eye



Model	Material		Bend Angle	Blade Thickness	Handle Length	Weight (Unit)	Secondary Package
	Blade	Handle					
MC-L-245	Stainless Steel (30Cr13)	ABS	Straight 25G	0.55 mm	117 mm	3.071 g	6 Pcs. / Box
MC-L-231	Stainless Steel (30Cr13)	ABS	Straight 20G	0.9 mm	117 mm	3.071 g	6 Pcs. / Box
MC-L-246	Stainless Steel (30Cr13)	ABS	Straight 25G	0.5 mm	117 mm	3.071 g	6 Pcs. / Box
MC-L-247	Stainless Steel (30Cr13)	ABS	Straight 23G	0.6 mm	117 mm	3.071 g	6 Pcs. / Box
MC-L-248	Stainless Steel (30Cr13)	ABS	Straight 21G	0.8 mm	117 mm	3.071 g	6 Pcs. / Box
MC-L-249	Stainless Steel (30Cr13)	ABS	Straight 19G	1.0 mm	117 mm	3.071 g	6 Pcs. / Box
MC-L-250	Stainless Steel (30Cr13)	ABS	Straight 18G	1.2 mm	117 mm	3.071 g	6 Pcs. / Box
MC-L-251	Stainless Steel (30Cr13)	ABS	Straight 16G	1.5 mm	117 mm	3.071 g	6 Pcs. / Box
MC-L-252	Stainless Steel (30Cr13)	ABS	45° 27G	0.4 mm	117 mm	3.071 g	6 Pcs. / Box
MC-L-253	Stainless Steel (30Cr13)	ABS	45° 26G	0.45 mm	117 mm	3.071 g	6 Pcs. / Box
MC-L-254	Stainless Steel (30Cr13)	ABS	45° 22G	0.7 mm	117 mm	3.071 g	6 Pcs. / Box
MC-L-255	Stainless Steel (30Cr13)	ABS	45° 21G	0.8 mm	117 mm	3.071 g	6 Pcs. / Box
MC-L-256	Stainless Steel (30Cr13)	ABS	45° 20G	0.9 mm	117 mm	3.071 g	6 Pcs. / Box
MC-L-257	Stainless Steel (30Cr13)	ABS	45° 19G	1.7 mm	117 mm	3.071 g	6 Pcs. / Box
MC-L-258	Stainless Steel (30Cr13)	ABS	45° 18G	1.8 mm	117 mm	3.071 g	6 Pcs. / Box
MC-L-259	Stainless Steel (30Cr13)	ABS	45° 16G	1.5 mm	117 mm	3.071 g	6 Pcs. / Box

Implant Slit Knife

Usage: to enlarge incision in cataract surgery



Model	Material		Bend Angle	Edge Width	Blade Thickness	Handle Length	Weight (Unit)	Secondary Package
	Blade	Handle						
KQ-L-278	Stainless Steel (30Cr13)	ABS	Implant Slit 45°, Bevel Up	2.8 mm	0.15 mm	117 mm	3.112 g	6 Pcs. / Box
KQ-L-280	Stainless Steel (30Cr13)	ABS	Implant Slit 45°, Bevel Up	3.0 mm	0.15 mm	117 mm	3.112 g	6 Pcs. / Box
KQ-L-281	Stainless Steel (30Cr13)	ABS	Implant Slit 45°, Bevel Up	3.5 mm	0.15 mm	117 mm	3.112 g	6 Pcs. / Box
KQ-L-282	Stainless Steel (30Cr13)	ABS	Implant Slit 45°, Bevel Up	4.1 mm	0.15 mm	117 mm	3.112 g	6 Pcs. / Box
KQ-L-283	Stainless Steel (30Cr13)	ABS	Implant Slit 45°, Bevel Up	5.1 mm	0.15 mm	117 mm	3.112 g	6 Pcs. / Box
KQ-L-284	Stainless Steel (30Cr13)	ABS	Implant Slit 45°, Bevel Up	5.2 mm	0.15 mm	117 mm	3.112 g	6 Pcs. / Box
KQ-L-284	Stainless Steel (30Cr13)	ABS	Implant Slit 45°, Bevel Up	5.6 mm	0.15 mm	117 mm	3.112 g	6 Pcs. / Box

Special Depth-limit Knife

Usage: corneal limbus incision while ensuring accurate cutting depth



Model	Material		Bend Angle	Edge Width	Blade Thickness	Weight (Unit)	Secondary Package
	Blade	Handle					
QG-L-260	Stainless Steel (30Cr13)	ABS	0°	1.3 mm	0.15 mm	3.071 g	6 Pcs. / Box
QG-L-261	Stainless Steel (30Cr13)	ABS	0°	1.4 mm	0.15 mm	3.071 g	6 Pcs. / Box
QG-L-262	Stainless Steel (30Cr13)	ABS	0°	1.5 mm	0.15 mm	3.071 g	6 Pcs. / Box
QG-L-263	Stainless Steel (30Cr13)	ABS	0°	1.6 mm	0.15 mm	3.071 g	6 Pcs. / Box
QG-L-264	Stainless Steel (30Cr13)	ABS	0°	1.7 mm	0.15 mm	3.071 g	6 Pcs. / Box
QG-L-265	Stainless Steel (30Cr13)	ABS	0°	1.8 mm	0.15 mm	3.071 g	6 Pcs. / Box
QG-L-266	Stainless Steel (30Cr13)	ABS	0°	1.9 mm	0.15 mm	3.071 g	6 Pcs. / Box
QG-L-267	Stainless Steel (30Cr13)	ABS	0°	2.0 mm	0.15 mm	3.071 g	6 Pcs. / Box
QG-L-268	Stainless Steel (30Cr13)	ABS	0°	2.2 mm	0.15 mm	3.071 g	6 Pcs. / Box
QG-L-269	Stainless Steel (30Cr13)	ABS	0°	2.3 mm	0.15 mm	3.071 g	6 Pcs. / Box
QG-L-270	Stainless Steel (30Cr13)	ABS	0°	2.4 mm	0.15 mm	3.071 g	6 Pcs. / Box
QG-L-271	Stainless Steel (30Cr13)	ABS	0°	2.5 mm	0.15 mm	3.071 g	6 Pcs. / Box
QG-L-272	Stainless Steel (30Cr13)	ABS	0°	2.65 mm	0.15 mm	3.071 g	6 Pcs. / Box
QG-L-273	Stainless Steel (30Cr13)	ABS	0°	2.75 mm	0.15 mm	3.071 g	6 Pcs. / Box
QG-L-274	Stainless Steel (30Cr13)	ABS	0°	2.8 mm	0.15 mm	3.071 g	6 Pcs. / Box
QG-L-275	Stainless Steel (30Cr13)	ABS	0°	2.85 mm	0.15 mm	3.071 g	6 Pcs. / Box
QG-L-276	Stainless Steel (30Cr13)	ABS	0°	3.0 mm	0.15 mm	3.071 g	6 Pcs. / Box
QG-L-277	Stainless Steel (30Cr13)	ABS	0°	3.2 mm	0.15 mm	3.071 g	6 Pcs. / Box
QG-L-278	Stainless Steel (30Cr13)	ABS	0°	2.2 mm	0.3 mm	3.071 g	6 Pcs. / Box
QG-L-279	Stainless Steel (30Cr13)	ABS	0°	2.2 mm	0.35 mm	3.071 g	6 Pcs. / Box
QG-L-280	Stainless Steel (30Cr13)	ABS	0°	2.2 mm	0.5 mm	3.071 g	6 Pcs. / Box
QG-L-281	Stainless Steel (30Cr13)	ABS	0°	2.2 mm	0.55 mm	3.071 g	6 Pcs. / Box
QG-L-282	Stainless Steel (30Cr13)	ABS	0°	2.2 mm	0.6 mm	3.071 g	6 Pcs. / Box
QG-L-283	Stainless Steel (30Cr13)	ABS	0°	2.2 mm	0.65 mm	3.071 g	6 Pcs. / Box


Special Golf Knife

Usage: pterygium excision




Model	Material		Bend Angle	Blade Length	Blade Width	Blade Thickness	Handle Length	Weight (Unit)	Secondary Package
	Blade	Handle							
BL-B-288	Stainless Steel (30Cr13)	ABS	Straight Golf Knife	25 mm	2.7 mm	1.2 mm	117 mm	3.071 g	6 Pcs. / Box

Product Features



Dura Blade


This blade is specially designed for dura mater. Its curved and blunt guard protects the tissues underneath during incision of dura mater or spinal dura mater.




Blade Thickness: 0.25 mm



DT-Z-352



Blade Thickness: 0.25 mm




DT-W-353

Arachnoid Blade

This blade is designed for arachnoid blade. The 0.3mm depth-limit design on the edge of circular blade prevents damage to deeper tissue.


The 2.0 mm edge width reveals excellent view of the incision.



Depth Controller (one side)

Blade Thickness: 0.25 mm


ϕ : 2.6 mm



DT-Z-354

Straight Blade

The ultra sharp 25° straight beveled blade is designed with a pointed tip at the front end. It is suitable for puncturing and cutting of micro-incision in narrow parts.




Blade Thickness: 0.38 mm
Angle: 25°



DT-Z-356

Curved Blade

The ultra sharp curved blade is with unique arc design can be used to hook tissue for cutting.




Blade Thickness: 0.38 mm



DT-Z-357

L-Hook Blade

The unique L-hook shaped blade with a small sickle design can be used to hook tissue for sharp piercing, and sharp cutting.



Blade thickness: 0.38 mm



DT-Z-358

Curved Scalpel Handle

310 mm



DB-W-363

210 mm



DB-W-362

185 mm



DB-W-361


Puncture Type Handle

200 mm




DB-W-369

100 mm



DB-W-368

Straight Scalpel Handle





NEUROTRONICS

Medical Solutions for Surgical Intervention
Neuro & Orthopedic Surgery

DISTRIBUTORS OF:

SCALPELS

EMG TUBES

IONM DEVICES

IONM CONSUMABLES

IMPLANTS & INSTRUMENT SETS



+27 31 562 0006



www.neurotronics.co.za



info@neurotronics.co.za